



Gas Chromatography Mass Spectrometry (GC-MS) an overview

Graeme Smith

Specialist Healthcare Practitioner

Willink Biochemical Genetics Lab

Manchester Centre for Genomic Medicine

Email: graeme.smith@cmft.nhs.uk



GCMS


- Basic Principle:
- Sample Prep
- Inlet
- Separation
- Detection



GCMS



GCMS



Sample Prep

- Liquid/Liquid
- Solid Phase
- Derivatisation:
 - » Volatile
 - » Thermally Stable



Gas Chromatography

- Inlet: Vapourise (Split and Splitless Injections)
- Separation

Inert stationary phase or GC Column
(Variety of properties, depending on the application)

Inert Carrier Gas, eg. He

Separates Primarily according to BP




Mass Spectrometry

Detector

Ion Source:


Most Applications use Electronic Ionisation
(EI)




Collision energy causes fragmentation



Ion Source



Quadrupole



Analysis Type

Scan

Scans all ions and masses

In order to avoid overloading the detector, this requires a split flow upon injection.



Analysis Type (Cont.)

Selected Ion Monitoring (SIM)

Only allows ions of a specific mass to pass through the quadrupole

Typically runs with a splitless injection

Can run multiple time windows and masses in the same run.



Quantitation

- Calibrant of known concentration e.g. 1mM
- Add Internal Standard.
- This generates a response, which in turn generates a specific concentration i.e. 1mM
- The Internal Standard is then added to each sample



Quantitation (Cont.)

- Concentration of the species in patient sample generates a response when measured with the internal standard.



Quantiation (Cont)

- Calibration

Single level or Multipoint

- Internal Standard

Stable Isotope Dilution or Non natural occurring species.



Applications

- Scan

Organic Acids

- SIM

VLCFA

Phytanic Acid

7DHC

MMA


Orotic Acid




Organic Acid Examples




Good Injection Profile




Chromatogram + MS




Mass Spectrum




Library Search



Bad Injection Profile



Really Bad Injection



Troubleshooting

- Maintenance
 - Inlet Septa
 - Inlet Liner
 - Column
 - Blocked Syringe
 - Insufficient Solvent




Further Troubleshooting

- Reduced Sensitivity
 - Air Leak
 - Dirty Ion Source


How do we know about this:



MS Tune



Air/Water Check




Air/Water Check

597x Air and Water Check

Instrument: Instrument #1
Mon Dec 29 13:02:11 2014 C:\MSDCHEM\1\5975\

Scan: 10.00 - 75.00 Samples: 8 Thresh: 0 Step: 0.10
60 peaks Base: 69.10 Abundance: 3287552



Current Params used: atune.u


Relative abundances:

18/69 = 21.72
28/69 = 6.94
32/69 = 2.48
44/69 = 3.65
28/18 = 31.95


Water% (Counts=714176)
Nitrogen% (Counts=228160)
Oxygen% (Counts=81456)
Carbon Dioxide% (Counts=119872)
Nitrogen/Water%




Benzoate




Valproate Interference




Abnormal 1



Abnormal 2



Abnormal 3




SIM Methods

Using VLCFA as an example



SIM Monitoring




Troubleshooting

- Ensure Time Windows are Correct
- If not run the standard on a scan method
- Reset Time Windows



Integration



Thank You For your
Time

Any Questions

